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GOOD-USE BIO

Broad Base Bio

FastTesting™ Plant genomic DNA extraction Kit

Cat# DE011

Kit Components

Component	Amount	
	50T	100T
Buffer A	15ml	30ml
Buffer B	5ml	10ml
Buffer BC	10ml	20ml
Buffer C	6.25ml	12.5ml
Buffer D	10ml	20ml
Spin column	50	100
Handbook	1 copy	1 copy

Important notes

1. Please add 21ml absolute ethanol to Buffer A and mix thoroughly before the first use.
2. Please operate in clean condition; All reagents should not be opened in the air for a long time.
3. All reagents, when stored at 15-25°C properly, are stable for 180 days.

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Protocol

1. Transfer up to 100mg tissue powder to a 1.5ml microcentrifuge tube.
2. Add 300ul A Buffer and 100 ul B Buffer, vortex gently for 5-10 minutes. If the sample is difficult to be lysised, please extend the incubation time.
3. Centrifuge for 3-5 minutes at 12,000 xg.
4. Transfer 200uL supernatant to a new 1.5ml tube, add 200uL BC buffer, incubate for 2 minutes.
5. Transfer the mixture to Spin Column, then centrifuge at 10,000 x g for 30 second. Discard flow-through.
6. Add 500µl C Buffer to the spin column. Centrifuge at 10,000 xg for 30 seconds. Discard flow-through.
7. Drying the spin for 5-10 minutes in a new 1.5ml tube.
8. Add 200ul TE buffer, Incubate for 2 minute, then centrifuge for 1 minute at 10,000 xg. The buffer in the tube contains the DNA.
9. The purified DNA can be used directly for kinds of downstream molecular biological experiments. Store at -20°C if not used immediately.

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